



S.A.L.T. Airway

Supraglottic Airway Laryngopharyngeal Tube

Frequently Asked Questions

Is the S.A.L.T. Airway approved by the FDA?

Yes, the S.A.L.T. Airway is registered with the FDA as a Class I Medical Device

Where can I purchase the S.A.L.T. Airway?

The S.A.L.T. Airway can only be purchased from one of our authorized distribution partners.

Contact your local

Ecolab – Alternative Healthcare Market, Region Managers for a complete listing

www.saltairway.com

How much does the S.A.L.T. Airway Cost?

The S.A.L.T. Airway Suggested Retail Price is \$18.95 each

Does the S.A.L.T. Airway come in multiple sizes?

Currently the S.A.L.T. Airway only comes in one size. We are hard at work developing additional sizes including a pediatric version

How do you size the S.A.L.T. Airway for a patient?

Corner of the mouth (proximal end) to the opposite angle of the jaw (distal end)

Who is currently using the S.A.L.T. Airway?

Many Public safety organizations around the country are using the S.A.L.T. Airway. Including; Major Fire Departments, EMS Departments, US Navy, US Army, Disaster Medical Response Teams, Hospitals and Surgical Centers

Can we obtain a sample for our Medical Director to evaluate?

Yes, visit http://www.mdimicrotek.com/salt_request_info.htm and complete the request form. Remember to note Medical Directors name and facility

Can we obtain a sample for our Training / Education facility?

Yes, visit http://www.mdimicrotek.com/salt_request_info.htm and complete the request form. Remember to note the facility name and contact person. Additionally please indicate total number of students and at which levels they are trained

Can the S.A.L.T. Airway be used in confined spaces with blind insertion?

Yes, we have received many stories of successful intubations from clinicians who said that their patient was in a confined space or trapped in an awkward position

What clinical research data is available?

We are currently conducting several large data collection studies. These studies are pre-hospital based departments with Medical Direction oversight. (Each study utilizes Paramedics trained in the use of the S.A.L.T. Airway. This is on going study; we are adding device uses on a regular basis. We will be releasing specific details very soon. We can confirm that the results have been extremely positive.)

